

For research use only. Not for use in diagnostic procedures.

## Prepito RNA Kit

*RNA purification from tissue samples*

Product no. CMG-2039

### Kit Components

<b>Magnetic Beads</b>	<b>Elution Buffer 6</b>
<b>Lysis Buffer 1</b>	<b>DNase Buffer</b>
<b>Binding Buffer 2</b>	<b>DNase I</b>
<b>Wash Buffer 3</b>	<b>Deep Well Plates</b>
<b>Wash Buffer 4</b>	<b>0.75 mL Reaction Tubes</b>
<b>Wash Buffer 5</b>	<b>Disposable Tips</b>

**Completion time:** approx. 110 minutes without sample preparation

### Equipment and other material to be provided by the user

1.5 mL reaction tube, disposable gloves, pipette and pipette tips with aerosol barrier (ensure that all used material is RNase free).

### Storage Conditions and Safety Information

Expiry dates are stated on the box and on each single component of the kit. Do not use any component of the kit beyond the expiration date. All kit components can be stored at room temperature.

Any further questions?


**chemagen Technology** technical support: +49 (0) 2401 805-501 | support.chemagen@perkinelmer.com



For research use only. Not for use in diagnostic procedures.

## Positioning Procedure

See “**Protocol Steps**” for detailed information



DNase position: 536  $\mu$ L MilliQ + 60  $\mu$ L DNase buffer + 4  $\mu$ L DNase I (10 U/ $\mu$ L)

Sample position: 600  $\mu$ L lysate

Pos. 4 second row for Disposable Tips; **! not used in this protocol !**

Pos. 3 Disposable Tips

Pos. 2 0.75 mL reaction tubes with 150  $\mu$ L **Magnetic Beads**

Pos. 1 0.75 mL reaction tubes with 70 – 150  $\mu$ L **Elution Buffer**

## Before You Start

- Check all kit components for integrity. In case of damages contact your supplier.
- Connect the tubes according to their numbering to the respective counterparts at the **chemagic 8-Pack**. Remove the lids from the individual buffer bottles in the **chemagic 8-Pack** and pierce the septum with the spike placed at the end of each tube. Place the **chemagic 8-Pack** upside down on the reagent holder and use the manual priming function for a complete filling of the dispensing system.
- RNase free water is needed and is not provided in this kit. Be sure that you have a sufficient volume for reconstitution of the **DNase I**.
- Dissolve the **DNase I** in the appropriate volume of sterile filtered 50 % glycerol/RNase free water solution (see flask label; glycerol is not provided in the kit).

Any further questions?

**chemagen Technology** technical support: +49 (0) 2401 805-501 | support.chemagen@perkinelmer.com



For research use only. Not for use in diagnostic procedures.

## Protocol Steps

1. Dissolve sample material in **Lysis Buffer 1** to obtain a total lysate volume of 600  $\mu$ L.
2. Carefully transfer 600  $\mu$ L lysate to each well of the deep well plate (DWP, riplate SW) defined as sample well (“Sample position: 600  $\mu$ L lysate”); see section above “Positioning Procedure”.
3. Prefill each well of the DWP defined as DNase well (“DNase position: 536  $\mu$ L MilliQ + 60  $\mu$ L DNase buffer + 4  $\mu$ L DNaseI (10 U/ $\mu$ L)”) with the respective DNase mixture; see section above “Positioning Procedure”.
4. Switch on the **chemagic Prepito** and wait for the self test to finish.
5. Press [**Change Protocol**].
6. Press [**Tissue**] in the Select Protocol Group window.
7. Select the **Prepito RNA Kit** protocol by pressing [**Prepito RNA Kit**] and confirm by pressing [**OK**].
8. Confirm the protocol selection in the Select Protocol Group window by pressing [**OK**].
9. Enter the 4 digit access code [**2272**] for authorization and confirm by pressing [**Enter**].
10. Press [**Start Process**].
11. Read the protocol information in the appearing information screen and confirm by pressing [**Continue**].
12. Select the sample positions and confirm by pressing [**OK**].
13. Enter the kit barcode with the barcode scanner and confirm by pressing [**OK**].
14. For the registration of the samples and the storage tubes press [**Yes**] and follow the instructions on the touch screen panel to enter the according barcodes.
15. Prepare the **chemagic Tip & Tube Rack** with the required materials. Place one 0.75 mL reaction tube filled with 70 – 150  $\mu$ L **Elution Buffer** (position 1), one 0.75 mL reaction tube filled with 150  $\mu$ L **Magnetic Beads** (position 2) and one **Disposable Tip** (position 3) for each sample into the positions according to the sample positions.

**!** *Shake the vessel with the Magnetic Beads vigorously until all Magnetic Beads are completely resuspended. An incomplete resuspension of the Magnetic Beads can cause a decreased yield of extracted nucleic acids.*

Any further questions?

**chemagen Technology** technical support: +49 (0) 2401 805-501 | support.chemagen@perkinelmer.com





For research use only. Not for use in diagnostic procedures.

16. An information screen indicates the previously selected sample positions. Ensure that the sample positions in the DWP correspond to the selected positions. Place the DWP on its default position on the tracking system and press [**Continue**].
17. Place the **chemagic Tip & Tube Rack** on its default position on the tracking system. Check for accurate fit of the DWP and **chemagic Tip & Tube Rack** and lock both by closing the safety latch.
18. Close the front door and start the automated isolation process by pressing [**Start**] immediately.

Any further questions?

**chemagen Technology** technical support: +49 (0) 2401 805-501 | support.chemagen@perkinelmer.com





For research use only. Not for use in diagnostic procedures.

## General Remarks

The **Elution Buffer** included in this kit is 10 mM Tris-HCl pH 8.0.

The **Magnetic Bead** suspension should be mixed vigorously before dispensing, otherwise the suspension is not homogenous and the DNA yield could be low.

## UV Measurements

In some cases you may find traces of **Magnetic Beads** left in the eluate. Such particles will not interfere with PCR and most downstream applications but may increase the background in UV measurements. In such a case, prior to UV analysis, we recommend an additional separation step using a manual separator (e.g. **chemagic Stand 2x12**, art. No. CMG-300) in order to separate any traces of particles.

Any further questions?

**chemagen Technology** technical support: +49 (0) 2401 805-501 | support.chemagen@perkinelmer.com

